

Estimation of Deleterious Particles in Lubricating Grease

ASTM D 1404

Product group(s): Grease

User group(s): Grease, Lubricating Grease

Scope: This test method describes a procedure for the detection and estimation of deleterious particles in lubricating grease.



Estimation Test Equipment

to be composed of:

- Deleterious Particles Determination Apparatus
- Test Plate

Technical Data

Dimensions: 30 x 25 x 25 cm

Weight: 3 kg

Main Unit

17-1530	Deleterious Particles Determination Apparatus ASTM D 1404
	<u>Consisting of:</u> body, test plate holders, loading screw, calibrated spring with scale for applying test load and removable cap assembly with milled slot and handle for rotating test plates.
17-1531	Test Plate made of plastic <u>Note!</u> 2 required for each test.

Order Guideline

Minimum equipment: 1x 17-1530, 2x 17-1531
 Spares (approx. 1 year): 17-1531
 Additional requirements: