

## Penetration Accessories - Rods & Rams

**Product group(s):** Penetration

**User group(s):** Adhesives, Asphalt, Automotive, Ball & Roller Bearings, Chemicals, Confectionery, Constituent, Cosmetics, Fats, Motor Vehicles, Penetration, Pharmaceuticals, Plastics, Polishes, Power Generation, Regulatory Authorities, Road Research, Rubber, Waste, Wax



### Penetration Accessories - Rods & Rams

More details are mentioned in the "Order Number" section

### Main Unit

<b>18-0260</b>	<b>Test Kit "Rod - Hollow"</b> <u>Suitable for:</u> semi-liquid fats, greases & pastes, emulsions, honey, paints & varnishes, potting compounds, ceramic pastes  <u>Consisting of:</u> 1 hollow rod, Ø 10 mm, 90°, 7 g, aluminum 3 interchangeable weights ( 3, 13 and 23 g)
<b>18-0261</b>	<b>Test Rod, Ø 10 mm, 10 g, plexiglas</b>
<b>18-0231</b>	<b>Test Cylinder, Ø 10 mm, 15 g, stainless steel</b>
<b>18-0232</b>	<b>Test Cylinder, Ø 15 mm, 35 g, stainless steel</b>
<b>18-0233</b>	<b>Test Cylinder, Ø 11 mm, 10 g, aluminum</b>
<b>18-0241</b>	<b>Test Ram AIB, 67 g, aluminum</b> for measurements on bread
<b>18-0243</b>	<b>DEUROBACK Test Cap, 50 g, stainless steel</b> for measurements on bread and confectionary
<b>18-2355</b>	<b>Test Ball, Ø 17 mm, 27.5 g, stainless steel</b>
<b>18-0255</b>	<b>Test Ball, Ø 6.3 mm, 6.3 g, brass</b> with shaft 100 mm long
<b>18-0251</b>	<b>Test Ram, Ø 8 mm, 15 g, aluminum</b> 40 mm effective length, with plane end face
<b>18-1726</b>	<b>Test Ram Ø 10 mm, 21.8 g, aluminum</b> 100 mm overall length, with plane end face
<b>18-0270</b>	<b>Test Cutter, 130 g, stainless steel</b>

18-0271	Test Ram, Ø 3 mm, 4.3 g, stainless steel
18-0272	Test Ram, Ø 6.3 mm, 6.4 g, stainless steel
18-0247	Platelike Test Ram "Adhesion", 10 cm <sup>2</sup> (Ø 35.7 mm)
18-0248	Platelike Test Ram "Adhesion", 20 cm <sup>2</sup> (Ø 50.5 mm)

### Order Guideline

Minimum equipment: Rods & Rams (see above)  
Additional requirements: